

# HealthView Integration Manager

LUMEDX's Integration Manager gives network administrators comprehensive control over their EHR/HIS and medical device interfaces through a single Web-based Management Console. The console uses browser technology to display the status of each interface, as well as a complete record of files imported and exported. Integration Manager's data management tools produce a log of unresolved files so administrators can remedy errors and upload corrected files into the system right away. In addition, Integration Manager features the industry's most extensive range of direct interfaces to cardiovascular and hospital information systems, so heart centers can connect to best-of-breed devices.

Integration Manager assures data integrity, streamlined workflow, and increased ROI from your cardiovascular IT investment.



## Features:

- For each interface view:
  - ◆ Status (running, stopped, etc.)
  - ◆ Number, date, and time of transactions for each interface
  - ◆ Last action
  - ◆ Number of successful, pending, and rejected transactions
- Remotely start and stop individual interfaces
- Resolve pending and rejected transactions
- Delete, edit, or re-register an interface
- Maintain log and intermediate database files

## Benefits

- Management Console is a single point-of-contact for all of your medical devices; access console from anywhere on your hospital network
- Interface to an extensive range of best-of-breed medical devices and hospital information systems
- Consolidate interfaces on one server; eliminate the expense of multiple communication or gateway servers
- Integration Manager works behind the scenes to streamline clinical and administrative workflow